

Two-Photon Direct Frequency Comb Spectroscopy of Potassium

M. E. Rowan, J. H. Baron, S. Chen and J. E. Stalnaker, *Department of Physics and Astronomy, Oberlin College*


Introduction

We discuss an experiment that uses direct frequency comb spectroscopy to study two-photon transitions in potassium. Atomic potassium is excited through two-photon transitions by use of the output of a stabilized optical frequency comb. The light generated by the comb is split, counter-propagated and focused into a heated vapor cell that contains potassium atoms. The repetition rate of the frequency comb is scanned and the potassium atoms are excited through various two-photon transitions. Transitions are detected via the fluorescence of the decaying excited state by use of a photomultiplier tube. We compare the experimental spectra with calculations of the two-photon transition amplitudes.


Optical Frequency Comb

- A mode-locked laser produces a series of ~30 fs long pulses with a repetition rate of 1 GHz.

A Series of Pulses:



- Phase coherence of the pulses leads to interference and the generation of an optical frequency comb.
- Dispersion causes a phase shift of the carrier wave relative to the envelope, resulting in a shift of the comb structure – the offset frequency f_0 .
- The frequency of each mode is given by an integer mode number, n , and two radio frequencies:
 1. The repetition rate, f_{rep} – the separation between modes.
 2. The offset frequency, f_0 – the shift of the comb relative to zero frequency.



Energy Level Diagram for Potassium

